

TEST REPORT T71-6-1

AD


QUALITATIVE ANALYSIS OF PRIMERS, TRACERS, IGNITERS,  
INCENDIARIES, BOOSTERS, AND DELAY COMPOSITIONS ON  
A MICRO SCALE BY USE OF INFRARED SPECTROSCOPY

by

David E. Chasan  
George Norwitz

June 1971

AMC CODE 4931.ON.6350



This document has been approved for public release and sale;  
its distribution is unlimited.

NATIONAL TECHNICAL  
INFORMATION SERVICE



**DEPARTMENT OF THE ARMY**  
**FRANKFORD ARSENAL**  
**Philadelphia, Pa. 19137**

## DISPOSITION

Destroy this report when it is no longer needed. Do not return it to the originator.

ADDITIONAL FOR		
CPSTI	WRITE SECTION	<input checked="" type="checkbox"/>
DOC	DIFF SECTION	<input type="checkbox"/>
UNANNOUNCED		<input type="checkbox"/>
JUSTIFICATION		
27		
DISTRIBUTION/AVAILABILITY CODES		
DMT	AYAN. AND/OR	SPECIAL
A		

## DISCLAIMER

The findings in this report are not to be construed as an official Department of the Army position, unless so designated by other authorized documents.

Security Classification

## DOCUMENT CONTROL DATA - R &amp; D

(Security classification of title, body of abstract and indexing annotation must be entered when the overall report is classified)

1. ORIGINATING ACTIVITY (Corporate author) Commanding Officer Frankford Arsenal ATTN: L3100-64-3, Philadelphia, Pa.		2a. REPORT SECURITY CLASSIFICATION UNCLASSIFIED	
		2b. GROUP	
3. REPORT TITLE Qualitative Analysis of Primers, Tracers, Igniters, Incendiaries, Boosters, and Delay Compositions on a Micro Scale by Use of Infrared Spectroscopy			
4. DESCRIPTIVE NOTES (Type of report and inclusive dates) Test Report			
5. AUTHOR(S) (First name, middle initial, last name) David E. Chasan George Norwitz			
6. REPORT DATE June 1971		7a. TOTAL NO. OF PAGES 35	7b. NO. OF REFS 3
8a. CONTRACT OR GRANT NO. AMC Code 4931.OM.6350		8b. ORIGINATOR'S REPORT NUMBER(S) T71- -1	
9. PROJECT NO. c. d.		9b. OTHER REPORT NO(S) (Any other numbers that may be assigned this report)	
10. DISTRIBUTION STATEMENT Release to OTS is authorized. This document has been approved for public release and sale; its distribution is unlimited.			
11. SUPPLEMENTARY NOTES		12. SPONSORING MILITARY ACTIVITY Army Materials Agency	
13. ABSTRACT The application of infrared spectroscopy to the detection of the constituents of primers, tracers, igniters, incendiaries, boosters, and delay compositions on a micro scale was investigated. It is shown that these constituents can be identified quickly and with certainty, using infrared pellet technique to detect organic and inorganic compounds and emission spectroscopy to identify the metals. In making the pellet, 1 to 2 mg of the material is ground with 300 mg of potassium bromide and the pellet formed in the die press. The grinding and pressing operation has been found to be completely safe even with the most sensitive explosives, once the material has been mixed with the potassium bromide. The infrared spectra of 43 of the most common ingredients of primers, tracers, igniters, incendiaries, boosters, and delay compositions are given over the range 2.5 to 50 microns. The qualitative analysis of seven typical compositions is demonstrated.			

DD FORM 1473

REPLACES DD FORM 1473, 1 JAN 64, WHICH IS OBSOLETE FOR ARMY USE.

Security Classification

14. KEY WORDS	LINK A		LINK B		LINK C	
	ROLE	WT	ROLE	WT	ROLE	WT
Primers Tracers Igniters Incendiaries Boosters Delay Compositions Infrared Spectroscopy						

T71-6-1

QUALITATIVE ANALYSIS OF PRIMERS, TRACERS, IGNITERS,  
INCENDIARIES, BOOSTERS, AND DELAY COMPOSITIONS ON  
A MICRO SCALE BY USE OF INFRARED SPECTROSCOPY

by

David E. Chasan  
George Norwitz

Project - Materials Testing Technology

AMS Code 4931.OM.6350  
PRON A1-0-P6350-01-AW-F6

This document has been approved for public release and sale;  
its distribution is unlimited.

Quality Assurance Directorate  
FRANKFORD ARSENAL  
Philadelphia, Pa. 19137

June 1971

## TABLE OF CONTENTS

	Page
I. SUMMARY . . . . .	1
II. RECOMMENDATIONS . . . . .	1
III. INTRODUCTION . . . . .	1
IV. STUDY	
A. Recommended Method for Preparation of Spectra . . .	2
B. Results and Discussion . . . . .	3
V. BIBLIOGRAPHY . . . . .	5
VI. TABLES . . . . .	6

## I. SUMMARY

The application of infrared spectroscopy to the detection of the constituents of primers, tracers, igniters, incendiaries, boosters, and delay compositions on a micro scale was investigated. It is shown that these constituents can be identified quickly and with certainty, using infrared pellet technique to detect organic and inorganic compounds and emission spectroscopy to identify the metals. In making the pellet, 1 to 2 mg of the material is ground with 300 mg of potassium bromide and the pellet formed in the die press. The grinding and pressing operation has been found to be completely safe even with the most sensitive explosives, once the material has been mixed with the potassium bromide. The infrared spectra of 43 of the most common ingredients of primers, tracers, igniters, incendiaries, boosters, and delay compositions are given over the range 2.5 to 50 microns. The qualitative analysis of seven typical compositions is demonstrated.

## II. RECOMMENDATIONS

It is recommended that the methods and spectra described in this report be included in MIL-STD-650 (Military Standard, Explosive: Sampling, Inspection and Testing).

## III. INTRODUCTION

The constituents of primers, tracers, igniters, incendiaries, boosters, and delay compositions are ordinarily identified by wet chemical tests for the organic compounds, anions, and metals (frequently emission spectrography is used for the metals). These wet tests are often untrustworthy, especially if applied on a micro scale.

This laboratory undertook an investigation on the application of infrared spectroscopy to the detection of the organic compounds and inorganic compounds present in primers, tracers, igniters, incendiaries, boosters, and delay compositions. It was believed that such an infrared technique together with emission spectroscopy (for inorganic materials that do not give infrared spectra) would provide an excellent means for the qualitative analysis of the constituents in question.

Pristera and Fredericks (3) have compiled spectra of most constituents of propellants and explosives. However, their compendium does not include many of the very sensitive explosive compounds found in primers, tracers, igniters, incendiaries, boosters, and delay compositions and it does not include the far infrared. Miller et al (1, 2) have compiled spectra of inorganic compounds but their collection does not include many of the inorganic compounds found in primers, tracers, igniters, incendiaries, boosters, and delay compositions. Pristera and Fredericks (3) and Miller et al (1, 2) used a prism type infrared spectrophotometer; the present authors used a grating type.

#### IV. STUDY

##### A. RECOMMENDED METHOD FOR PREPARATION OF SPECTRA

###### 1. Apparatus and Reagents

Perkin-Elmer Model 621 Infrared Spectrophotometer.

Evacuatable potassium bromide die (Model 186-0025, Perkin-Elmer Corp., Norwalk, Conn).

Adapters for mounting KBr discs on spectrophotometer.

Twenty-ton hydraulic press.

Vacuum pump.

Mortars (o.d. 50 mm) with pestles.

Plastic microspatula.

Potassium bromide, infrared grade (Harshaw Chemical Co., Cleveland, Ohio).



## 2. Preparation of Potassium Bromide Pellets

Grind 1 - 2 grams of infrared grade KBr (in small portions) and store in a desiccator until needed.

Add 1 - 2 mg of sample to the mortar with a plastic microspatula (this is most conveniently done by first taring the mortar). Weigh 300 mg of the preground KBr into an aluminum balance pan and add all the KBr to the mortar. Using the plastic microspatula, thoroughly mix the KBr with the sample. Grind the sample-KBr mix for about three minutes. Scrape the sample-KBr mix loose from the mortar surface with the plastic microspatula and transfer it to the die. Level the sample-KBr mix with the microspatula and insert the plunger. Evacuate the die for three minutes with a vacuum pump without applying pressure, then apply 18,000 pounds pressure with a hydraulic press for five minutes.

Disassemble the die and push out the KBr pellet with a gentle application of the hydraulic press. Mount the pellet on an adapter and scan the spectrum from 2.5 to 50 microns with a blank KBr pellet in the reference beam. Use a dry air purge for the far infrared region (30 to 50 microns). The following settings were used on the Perkin-Elmer Model 621 infrared spectrophotometer: slit program, 1000; gain 5.1; attenuator speed, 1100; scan drive, 1; scan time, 16; suppression, 8; scale, 1x; source current, 0.8.

## B. RESULTS AND DISCUSSION

The spectra obtained for the 43 commonly used ingredients of primers, tracers, igniters, incendiaries, boosters, and delay compositions are shown in Figures 1 to 43. These ingredients are arranged (Table I) according to function as follows: primary explosives, high explosives, oxidizers, color intensifiers, binders, and miscellaneous. Compounds having more than one usage are listed according to their major function; for example, polyvinyl chloride is used as a binder or color intensifier but is listed under binders. Lead thiocyanate, the only fuel that gives a spectrum, is listed under miscellaneous; other fuels (magnesium, aluminum, zirconium, titanium, boron, silicon, and metallic hydrides) do not give infrared spectra.

The spectra of seven typical explosive compositions together with the notation of the peaks identifying the principle constituents are shown in Figures 44 to 50. The conclusions concerning the analysis of these compositions drawn from these spectra and the emission spectrographic analysis are shown in Table II.

It is seen that normal lead styphnate can be distinguished from basic lead styphnate by infrared. The different oxides of lead, iron, and barium may also be distinguished ( $\text{PbO}$ ,  $\text{Pb}_3\text{O}_4$ ,  $\text{Fe}_2\text{O}_3$ ,  $\text{BaO}$ , and  $\text{BaO}_2$  give unique spectra, while  $\text{PbO}$  and  $\text{Fe}_3\text{O}$  give no spectra).

In some spectra there is evidence of light scattering in the high frequency end of the spectrum. This is caused by inadequate grinding or the presence of such ingredients as antimony sulfide, magnesium, aluminum or zirconium. However, in no instance was the problem severe enough to interfere with the qualitative use of the spectrum.

The majority of spectra indicated the presence of small amounts of moisture; hence, care must be taken to avoid misinterpreting peaks occurring in the regions of  $3420\text{ cm}^{-1}$  and  $1630\text{ cm}^{-1}$ .

A special concentration procedure must be used to detect tetracene (2 to 4%) present in some lead styphnate primers. In this concentration procedure, about 50 mg of the primer contained in a sintered glass crucible is washed successively with ammonium acetate solution (20%), water, and acetone to remove lead styphnate, barium nitrate, and PETN. About 3 or 4 mg of the residue, which contains the tetracene and antimony sulfide is then used for the pellet.

It might be thought that preparing potassium bromide pellets from very sensitive explosive materials would be hazardous. This did not prove to be the case for the following reasons: (a) a very small amount of the explosive is used, (b) the potassium bromide acts as a diluent.

In the emission spectrographic analysis of explosives, it is customary in this laboratory to first destroy the explosive by treatment with concentrated nitric acid and evaporation to dryness.

## V. BIBLIOGRAPHY

- (1) Miller, F. A., Carlson, G. L., Bently, F. F., and Jones, W. H., Spectrochim. Acta, 15, 135 (1960).
- (2) Miller, F. A., and Wilkens, C. H., Anal. Chem., 24, 1253 (1952).
- (3) Pristera, F., and Fredricks, W., Picatinny Arsenal Technical Memorandum 1887, September 1969.

## VI. TABLES

Table I. Identification of Spectra

<u>Spectrum</u>	<u>Code</u>
1. Lead azide	Primary Explosive 1
2. Mercury fulminate	Primary Explosive 2
3. Tetracene	Primary Explosive 3
4. Normal lead styphnate	Primary Explosive 4
5. Basic lead styphnate	Primary Explosive 5
6. RDX (cycotrimethylenetrinitramine)	High Explosive 1
7. HMX (cyclotetramethylenetetranitramine)	High Explosive 2
8. Tetryl (2,4,6-trinitrophenylmethylnitramine)	High Explosive 3
9. PETN (pentaerythrite tetranitrate)	High Explosive 4
10. TNT (2,4,6-trinitrotoluene)	High Explosive 5
11. Nitrocellulose	High Explosive 6
12. Styphnic acid	High Explosive 7
13. Ammonium nitrate	Oxidizer 1
14. Sodium nitrate	Oxidizer 2
15. Potassium nitrate	Oxidizer 3
16. Barium nitrate	Oxidizer 4
17. Strontium nitrate	Oxidizer 5
18. Lead nitrate	Oxidizer 6
19. Potassium chlorate	Oxidizer 7
20. Ammonium perchlorate	Oxidizer 8

Table I. Identification of Spectra (cont'd)

<u>Spectrum</u>	<u>Code</u>
21. Potassium perchlorate	Oxidizer 9
22. Barium chromate	Oxidizer 10
23. Lead chromate	Oxidizer 11
24. Barium peroxide	Oxidizer 12
25. Strontium peroxide	Oxidizer 13
26. Lead monoxide	Oxidizer 14
27. Red lead ( $Pb_3O_4$ )	Oxidizer 15
28. Ferric oxide	Oxidizer 16
29. Molybdenum trioxide	Oxidizer 17
30. Dechlorane (perchloropentacyclodecane)	Color Intensifier 1
31. Sodium oxalate	Color Intensifier 2
32. Strontium oxalate	Color Intensifier 3
33. Gum arabic	Binder 1
34. Polyvinyl chloride	Binder 2
35. Parlon (chlorinated rubber)	Binder 3
36. Stearic acid	Binder 4
37. Wax	Binder 5
38. Calcium resinate	Binder 6
39. Ethyl cellulose	Binder 7
40. Oxamide	Misc. 1

Table I. Identification of Spectra (cont'd)

<u>Spectrum</u>	<u>Code</u>
41. Barium monoxide	Misc. 2
42. Lead thiocyanate	Misc. 3
43. Calcium carbonate	Misc. 4
44. Detonator	Mix 1
45. Igniter	Mix 2
46. Primer	Mix 3
47. Primer	Mix 4
48. Primer	Mix 5
49. Detonator	Mix 6
50. Detonator	Mix 7

Table II. Analysis of Typical Compositions.

<u>Spectrum</u>	<u>Infrared</u>	<u>Emission Spectroscopy</u>	<u>Combined Analysis</u>
44	RDX TNT		RDX TNT
45	BaCrO <sub>4</sub> Pb <sub>3</sub> O <sub>4</sub>	Pb, Cr, Ba, Zr	BaCrO <sub>4</sub> Pb <sub>3</sub> O <sub>4</sub> Zr
46	Hg Fulminate KClO <sub>3</sub>	K, Sb, Hg	Hg Fulminate KClO <sub>3</sub> Sb <sub>2</sub> S <sub>3</sub>
47	Hg Fulminate KClO <sub>3</sub> Ba(NO <sub>3</sub> ) <sub>2</sub>	K, Ba, Sb, Hg	Hg Fulminate KClO <sub>3</sub> Ba(NO <sub>3</sub> ) <sub>2</sub> Sb <sub>2</sub> S <sub>3</sub>
48	Normal Pb Styphnate PETN Ba(NO <sub>3</sub> ) <sub>2</sub>	Al, Sb, Pb, Ba	Normal Pb Styphnate PETN Ba(NO <sub>3</sub> ) <sub>2</sub> Al Sb <sub>2</sub> S <sub>3</sub>
49	Pb Azide PETN	Pb	Pb Azide PETN
50	RDX Stearic Acid		RDX Stearic Acid

